

LeXBot Flex


Fully-automated NGS Workstation

Introduction


LeXBot Flex Fully-automated NGS Workstation, independently developed by LexigenBio, is designed for automated library preparation in gene sequencing. It supports a wide range of automated operations, including sample dispensing, nucleic acid extraction, DNA normalization, and sample mixing, library preparation, and hybrid capture, significantly improving laboratory efficiency. The workstation features four independent pipetting channels and seamlessly integrates a plate-gripping module, thermal cycler (brand optional), temperature control module, heating and oscillating module, purification module, fluorometer module, UV lamp set, and HEPA filtration system. With 20 built-in SBS standard plate positions, this workstation achieves a maximum throughput of up to 32 library preparations or 16 hybrid captures per run.



Application



Feature



LeXBot Flex


Flexible Integration

- 20 standard SBS built-in plate positions, allowing modules to be freely configured to meet diverse experimental needs.
- Seamless integration of multifunctional components, including thermal cycler, temperature control, heating and oscillation, purification, and fluorometer modules.
- A 3D robotic arm with a plate gripper intelligently handles various consumables, ensuring complete deck coverage.

Fully-automated NGS Workstation

versatile Application

- Four independent pipetting channels with liquid level detection and intelligent tracking for precise liquid handling.
- Each channel supports variable spacing, volume, and height pipetting to accommodate diverse requirements.
- Flexibly supports 1-32 library preparations and 1-16 hybridization captures, meeting varying throughput demands.



LeXBot Flex


Fully-automated NGS Workstation

Secure and Reliable

- Bright, multi-color status indicator lights provide real-time updates on operational status, minimizing downtime caused by hardware or software errors.
- Equipped with an emergency pause button for one-touch emergency stops, ensuring human-machine safety.
- Standard HEPA filtration system and UV lamp maintain a continuously clean operating environment.



Dynamic Operating Deck



Functional Component



Temperature Control Module

- Consists of 2 temperature control modules.
- Temperature control range 4°C to 75°C.
- Temperature accuracy $\leq \pm 1.0^\circ\text{C}$.
- Temperature uniformity $\leq 1.0^\circ\text{C}$.
- Temperature fluctuation $\leq \pm 0.5^\circ\text{C}$.



Heating and Oscillating Module


- Temperature control range: (room temperature +5°C) to 99°C.
- Temperature accuracy $\leq \pm 2.0^\circ\text{C}$.
- Temperature uniformity $\leq 1.5^\circ\text{C}$.
- Temperature fluctuation $\leq \pm 0.5^\circ\text{C}$.
- Oscillating speed adjustable from 200 rpm to 2000 rpm.
- Oscillating speed accuracy $\leq \pm 25$ rpm.



Thermal Cycler

- Brand optional.
- Temperature control range: 4°C to 99°C.
- Heated lid temperature range: 30°C to 110°C.
- Average heating rate $\geq 3.0^\circ\text{C}/\text{sec}$.
- Average cooling rate $\geq 1.5^\circ\text{C}/\text{sec}$.
- Temperature accuracy $\leq \pm 0.5^\circ\text{C}$.
- Temperature uniformity $\leq 1.0^\circ\text{C}$.
- Temperature fluctuation $\leq \pm 0.2^\circ\text{C}$.

Purification Module




- Includes one liftable 96-well deep-well plate ring magnetic stand and one movable 96-well PCR plate ring magnetic stand.



Fluorometer Module


- Simultaneous concentration measurement for up to 8 samples in one movement, with a maximum of 32 samples per run.
- Capable of performing automatic normalization and pooling operations.
- Fast: accurate concentration measurement within 3-5 sec.
- Ultra-low input requirement: only 1-20 μL per sample.



PCR Plate Carrier



Composite Plate Carrier



Reagent Reservoir Rack



Deep Well Plate Carrier




Four Independent Pipetting Channels

- Supports variable spacing, volume, and height for pipetting.
- Equipped with liquid-level detection, enabling intelligent liquid-level tracking and segmented Z-axis descent (movement, detection, aspiration).
- Pipetting range: 1 μL -1000 μL .
- Performance specifications:
 1 μL - CV < 8%, accuracy < $\pm 10\%$.
 200 μL - CV < 1%, accuracy < $\pm 2\%$.
 1000 μL - CV < 1%, accuracy < $\pm 1\%$.


Consumables



**1000 μL
filtered tip box**



**200 μL
filtered tip box**



**50 μL
filtered tip box**



**1.3 mL U-bottom
deep well plate**



**fully skirted
96-well PCR plate**




**25 mL
reagent reservoir**



**0.5 mL/2 mL
cryogenic tube**



**0.2 mL thin-wall
PCR 8-strip tube**



**PCR plate
sealing pad**



**fluorometer
8-strip tube**

Specification

Host parameters	Size	1160 mm (length) × 815 mm (width) × 970 mm (height)		weight	180 KG
	Maximum output power	1700 VA			
power requirements	Voltage	AC 100-240V		Frequency	50/60 Hz
Operating environment	Temperature range	10°C ~ 30 °C		Relative humidity 20% ~ 80% , non-condensing	
	Barometric pressure range	86 kPa ~ 106 kPa			
Robot arm	Repeated positioning	±0.1 mm			

Ordering Information

catalog	product	Scale
LX_LeXBotFlex	LeXBot Flex Fully-automated NGS Workstation	RU0

Statement

For research use only. Not for use in diagnostic procedures.

Without the written permission of LexigenBio, no other individual or organization may reproduce, copy, edit, expand, reduce or translate the contents of this document in any form into other languages for any purpose.

In case of improper use, LexigenBio Inc reserves all rights