

Hematologic Malignancy Detection | LeXHema Panel v2.0

Background

Hematologic malignancies are a group of clonal malignant diseases originating from the hematopoietic system, characterized by high heterogeneity. They include major types such as leukemia, lymphoma, multiple myeloma, and myelodysplastic syndromes. Different hematologic malignancies exhibit significant variations in their molecular genetic features, involving a range of genetic alterations such as base substitutions, insertions/deletions, copy number variations (CNVs), and gene fusions. Consequently, comprehensive and precise molecular detection is crucial for the accurate diagnosis, targeted treatment strategy formulation, and minimal residual disease (MRD) monitoring of these disorders.

Traditional detection techniques—such as FISH, qPCR, and Sanger sequencing—offer certain advantages for detecting specific mutations. However, they often cover only a limited set of genes or known hotspot mutations, making it challenging to fully decipher the complex spectrum of genetic abnormalities. This is particularly problematic for detecting low-frequency mutations and discovering novel fusion genes. The emergence of high-throughput sequencing (NGS) technology has delivered a breakthrough, enabling panoramic molecular profiling of hematologic malignancies.

DNA-based detection can cover coding regions as well as portions of intronic regions, allowing for comprehensive analysis of point mutations, indels, CNVs, and fusion breakpoints. This approach is well-suited for screening both known and potential pathogenic variants. However, gene fusion detection via DNA methods relies on intronic probes, which can be compromised by the extensive length and repetitive nature of intronic regions, potentially resulting in missed fusion events. In contrast, RNA-based detection directly targets fusion gene detection at the transcript level. It not only precisely identifies genuine fusion events and their partner genes but also evaluates whether these fusions result in aberrant transcripts, thereby elucidating their functional impact. Thus, the integration of both DNA and RNA workflows achieves a multidimensional analysis of genetic alterations: DNA detection provides broad coverage for mutation screening, while RNA detection accurately confirms fusion expression and functionality. Together, they ensure comprehensive detection and enhanced accuracy of results.

Introduction

LeXHema Panel v2.0 targets 481 genes associated with hematologic malignancies, comprehensively covering a wide range of genetic alterations including base substitutions, insertions/deletions, amplifications, and gene fusions.

This panel consists of two components:

LeXHema Panel v2.0_DNA Tube: Covers 1.7 Mb of genomic regions, including full coding regions of 436 genes, partial intronic regions of 12 genes, and the *IGH*Switch regions.

LeXHema Panel v2.0_RNA Tube: Covers 146 genes, optimized for fusion gene detection.

By integrating multi-omics approaches, employing a highly targeted design, and enhancing clinical utility, **LeXHema Panel v2.0** offers an efficient and cost-effective solution for precision diagnosis, treatment planning, and longitudinal management of hematologic malignancies.

Feature

Comprehensive Multi-type Mutation Coverage

- Encompasses 481 hematologic malignancy-related genes, enabling thorough analysis of base substitutions, indels, amplifications, and gene fusions.

Dual Workflow for Enhanced Accuracy

- RNA workflow compensates for the limitations of the DNA workflow by accurately identifying fusion gene partners and validating the functional impact of fusion events. The synergistic use of both workflows ensures comprehensive and precise detection.

Flexible Customization

- Supports flexible subdivision into multiple sub-panels or the addition of target genes, allowing adaptation to diverse hematologic malignancy diagnostic needs in both research and clinical applications.

Performance

LeXHema Panel v2.0_DNA

Capture Performance




Figure 1. Capture performance of LeXHema Panel v2.0_DNA applied to various types of gDNA samples. **A.** Mappability & On-target rate; **B.** Target covered; **C.** Average sequencing depth after deduplication; **D.** GC bias. Using 50 ng of various gDNA samples, pre-libraries were prepared with the LexPrep EZ DNA Library Preparation Kit coupled with the LexPrep Universal Stubby Adapter (UDI) Module, followed by hybrid capture using LeXHema Panel v2.0_DNA with LexPrep Hybrid Capture Reagents. Sequencing was performed on both NovaSeq 6000 (PE150) and DNBSEQ-T7 (PE150) platforms, with a random subset of 1.3 Gb used for data analysis.

Note: LDT-831: hematologic malignancy gDNA quality control (LDT Bioscience; LDT-831A: AF = 5%; LDT-831D: AF = 0%). Male: human genomic DNA standard (Promega, G1471).

Stable Detection of Multiple Mutation Types

Table 1. Detection of various mutation types in different standard samples using LeXHema Panel v2.0_DNA.

Type	Variant	LDT-831A			LDT-831D			Male		
		Reference	Observed		Reference	Observed		Reference	Observed	
			NovaSeq	DNBSEQ		NovaSeq	DNBSEQ		NovaSeq	DNBSEQ
SNV	<i>JAK2</i> V617F	5%	4.4%	4.9%	0%	0%	0%	0%	0%	0%
SNV	<i>KIT</i> D816V	5%	5.1%	3.6%	0%	0%	0%	0%	0%	0%
SNV	<i>TP53</i> R175H	5%	3.5%	6.4%	0%	0%	0%	0%	0%	0%
SNV	<i>RUNX1</i> R201Q	5%	5.7%	4.0%	0%	0%	0%	0%	0%	0%
SNV	<i>IDH1</i> R132H	5%	5.9%	4.3%	0%	0%	0%	0%	0%	0%
Insertion	<i>FLT3</i> E598_Y599insFDFREYE	5%	3.8%	3.1%	0%	0%	0%	0%	0%	0%
Fusion	<i>BCR-ABL1</i> (e14a2)	5%	4.4%	3.4%	0%	0%	0%	0%	0%	0%
Fusion	<i>CCDC6-RET</i>	5%	12.5%	20.2%	0%	0%	0%	0%	0%	0%

LeXHema Panel v2.0_RNA

Capture Performance




Table 2. Detection of fusion genes in K562 cell line RNA using LeXHema Panel v2.0_RNA.

Fusion type	Platform	Split reads1	Split reads2	Discordant mates	Coverage1	Coverage2
<i>BCR-ABL1</i>	NovaSeq	298	214	14	948	882
	DNBSEQ	223	204	10	852	849
<i>NUP214-XKR3</i>	NovaSeq	86	28	6	1114	30
	DNBSEQ	87	23	6	1155	27

Note: Data are reported per M read pairs (approximately 0.3 Gb); analysis was performed using Arriba.

Gene List

LeXHema Panel v2.0_DNA_436 genes_Full CDS

ABL1	ABL2	ACTB	AKT1	AKT2	AKT3	ALK	AMER1	ANKRD26	AP3B1	APC	AR	ARAF	ARHGAP26	ARID1A
ARID1B	ARID2	ARID5B	ASXL1	ASXL2	ASXL3	ATG5	ATM	ATR	ATRX	ATXN1	AURKA	AURKB	AXIN1	AXL
B2M	BAALC	BAP1	BARD1	BCL10	BCL11A	BCL11B	BCL2	BCL6	BCL7A	BCOR	BCORL1	BIRC3	BLM	BLNK
BRAF	BRCA1	BRCA2	BRD4	BTG1	BTG2	BTX	BTLA	CALR	CARD11	CBFB	CBL	CBLB	CCND1	CCND2
CCND3	CCNE1	CCR4	CCR7	CCT6B	CD22	CD27	CD274	CD28	CD58	CD70	CD79A	CD79B	CD83	CDC42
CDC73	CDH1	CDK12	CDK4	CDK6	CDK8	CDKN1A	CDKN1B	CDKN1C	CDKN2A	CDKN2B	CDKN2C	CEBPA	CEP57	CHD2
CHD8	CHEK1	CHEK2	CIITA	CKS1B	CREBBP	CRLF2	CSF1R	CSF3R	CTCF	CTLA4	CTNNB1	CUX1	CXCR4	CYLD
DAXX	DDR2	DDX3X	DDX41	DHX15	DNM2	DNMT3A	DOT1L	DTX1	DUSP2	DUSP22	EBF1	ECT2L	EED	EGFR
EP300	EPCAM	EPHA2	EPHA3	EPHA5	EPHA7	EPOR	ERBB2	ERBB3	ERBB4	ERCC4	ERG	ESR1	ETNK1	ETS1
ETV1	ETV4	ETV6	EZH2	FAM46C	FANCA	FANCC	FANCD2	FANCE	FANCF	FANCG	FANCI	FANCL	FANCL	FANCM
FAS	FAT1	FAT4	FBXO11	FBXW7	FGFR1	FGFR2	FGFR3	FGFR4	FH	FLCN	FLT1	FLT3	FLT4	FOXO1
FOXO3	FOXP1	FYN	GADD45B	GATA1	GATA2	GATA3	GNAI1	GNAI3	GNAQ	GNAS	GNB1	GRIN2A	HACE1	HBA1
HBA2	HBB	HDAC1	HDAC4	HDAC7	HGF	HIST1H1A	HIST1H1C	HIST1H1D	HIST1H1E	HNF1A	HRAS	ID3	IDH1	IDH2
IFNAR2	IGF1R	IGH	IGK	IGL	IGLL5	IKBKE	IKZF1	IKZF2	IKZF3	IL3	IL7R	INPP4B	INPP5D	IRF1
IRF4	IRF8	ITK	ITPKB	JAK1	JAK2	JAK3	JARID2	JUN	KAT6A	KDM2B	KDM5A	KDM5C	KDM6A	KDR
KIT	KLF2	KLHL5	KLHL6	KMT2A	KMT2B	KMT2C	KMT2D	KRAS	LEF1	LMO2	LRP1B	LTB	LYN	LYST
MAF	MAFB	MALT1	MAP2K1	MAP2K2	MAP2K4	MAP3K1	MAP3K14	MAPK1	MCL1	MDM2	MDM4	MED12	MEF2B	MEF2D
MEN1	MET	MFHAS1	MGA	MITF	MLH1	MPL	MSH2	MSH3	MSH6	MTOR	MUTYH	MYB	MYC	MYCL
MYCN	MYD88	NBN	NCKAP1L	NCSTN	NF1	NF2	NFE2L2	NFKB1	NFKB2	NFKBIA	NFKBIE	NKX2-1	NLR4	NOTCH1
NOTCH2	NPM1	NRAS	NSD1	NTSC2	NTRK1	NTRK3	NUP98	P2RY8	PAG1	PAK3	PALB2	PAX5	PBRM1	PC
PDCD1	PDCD1LG2	PDGFRA	PDGFRB	PDK1	PHF6	PIK3CA	PIK3CD	PIK3R1	PIK3R2	PIM1	PLCG1	PLCG2	PML	PMS2
POLE	POT1	POU2AF1	PPM1D	PPP2R1A	PRDM1	PRF1	PRKCB	PRPF8	PTCH1	PTEN	PTPN1	PTPN11	PTPN2	PTPN6
PTPRD	PTPRK	PTPRO	RAB27A	RAD21	RAD50	RAD51	RAF1	RARA	RASGEF1A	RB1	RECQL4	REL	RELN	RET
RHOA	RICTOR	RNF43	ROS1	RPL10	RPS15	RPTOR	RUNX1	SAMD9	SAMD9L	SBD5	SDHA	SDHB	SDHC	SDHD
SETBP1	SETD2	SF3B1	SGK1	SH2B3	SH2D1A	SLC19A1	SLC7A7	SMAD2	SMAD4	SMARCA2	SMARCA4	SMARCB1	SMC1A	SMC3
SMO	SOC1	SOX11	SOX2	SPEN	SPOP	SRC	SRP72	SRSF2	STAG2	STAT3	STAT4	STAT5A	STAT5B	STAT6
STK11	STX11	STXB2	SUFU	SUZ12	SYK	TBL1XR1	TCF3	TCL1A	TERC	TERT	TET1	TET2	TET3	TGFB2
TNFAIP3	TNFRSF11A	TNFRSF14	TNFRSF17	TNFRSF18	TOP1	TOP2A	TP53	TP63	TP73	TRAF2	TRAF3	TRAF5	TSC1	TSC2
TSHR	TYK2	U2AF1	UNC13D	VAV1	VHL	WHSC1	WT1	XIAP	XPC	XPO1	ZAP70	ZBTB7A	ZEB2	ZMYM3
ZRSR2														

LeXHema Panel v2.0_DNA_12 genes_Selected Intron

ALK	BCL6	BCR	BRAF	ETV6	JAK2	KMT2A	MYC	PDGFRA	PDGFRB	RARA	RET
-----	------	-----	------	------	------	-------	-----	--------	--------	------	-----

(Continued)

LeXHema Panel v2.0_RNA_146 genes

<i>ABL1</i>	<i>ABL2</i>	<i>ALK</i>	<i>ARHGAP26</i>	<i>ARID1A</i>	<i>ASXL1</i>	<i>ATG5</i>	<i>BCL10</i>	<i>BCL11A</i>	<i>BCL11B</i>	<i>BCL2</i>	<i>BCL3</i>	<i>BCL6</i>	<i>BCL7A</i>	<i>BCL9</i>
<i>BCOR</i>	<i>BCR</i>	<i>BIRC3</i>	<i>BRAF</i>	<i>BTG1</i>	<i>CBFA2T3</i>	<i>CBFB</i>	<i>CBL</i>	<i>CCND1</i>	<i>CCND2</i>	<i>CCND3</i>	<i>CD274</i>	<i>CDK6</i>	<i>CEBPE</i>	<i>CIITA</i>
<i>CREBBP</i>	<i>CRLF2</i>	<i>CTNNB1</i>	<i>DEK</i>	<i>DUSP22</i>	<i>DUX4</i>	<i>EGFR</i>	<i>EP300</i>	<i>EPOR</i>	<i>ERBB2</i>	<i>ERG</i>	<i>ETSI</i>	<i>ETV1</i>	<i>ETV4</i>	<i>ETV6</i>
<i>EWSR1</i>	<i>FGFR1</i>	<i>FGFR2</i>	<i>FGFR3</i>	<i>FOXO1</i>	<i>FOXO3</i>	<i>FOXP1</i>	<i>GLIS2</i>	<i>HLF</i>	<i>IGH</i>	<i>IGK</i>	<i>IgL</i>	<i>IKZF1</i>	<i>IL3</i>	<i>IRF4</i>
<i>ITK</i>	<i>JAK1</i>	<i>JAK2</i>	<i>JAK3</i>	<i>KAT6A</i>	<i>KIF5B</i>	<i>KMT2A</i>	<i>LMO1</i>	<i>LMO2</i>	<i>LYN</i>	<i>MAF</i>	<i>MAFB</i>	<i>MALT1</i>	<i>MECOM</i>	<i>MEF2D</i>
<i>MLLT1</i>	<i>MLLT10</i>	<i>MLLT3</i>	<i>MIN1</i>	<i>MTCP1</i>	<i>MYB</i>	<i>MYC</i>	<i>MYH11</i>	<i>NCOA2</i>	<i>NCOR1</i>	<i>NF1</i>	<i>NF2</i>	<i>NFKB2</i>	<i>NOTCH1</i>	<i>NPM1</i>
<i>NSD1</i>	<i>NTRK1</i>	<i>NTRK3</i>	<i>NUP214</i>	<i>NUP38</i>	<i>NUTM1</i>	<i>P2RX8</i>	<i>PAX5</i>	<i>PBX1</i>	<i>PDCD1LG2</i>	<i>PDGFRA</i>	<i>PDGFRB</i>	<i>PICALM</i>	<i>PIM1</i>	<i>PML</i>
<i>POU2AF1</i>	<i>PRDM1</i>	<i>PTCH1</i>	<i>RAF1</i>	<i>RARA</i>	<i>RARB</i>	<i>RARG</i>	<i>RBM15</i>	<i>RET</i>	<i>ROS1</i>	<i>RUNX1</i>	<i>RUNX1T1</i>	<i>STAT6</i>	<i>STIL</i>	<i>SYK</i>
<i>TAF15</i>	<i>TAL1</i>	<i>TBL1XR1</i>	<i>TCF3</i>	<i>TCL1A</i>	<i>TCL1B</i>	<i>TET1</i>	<i>TLX1</i>	<i>TLX3</i>	<i>TNFRSF11A</i>	<i>TOP1</i>	<i>TP63</i>	<i>TPM3</i>	<i>TPM4</i>	<i>TRA</i>
<i>TRB</i>	<i>TRD</i>	<i>TRG</i>	<i>TYK2</i>	<i>VAV1</i>	<i>WHSC1</i>	<i>ZBTB16</i>	<i>ZEB2</i>	<i>ZNF292</i>	<i>ZNF362</i>	<i>ZNF384</i>				

Ordering Information

Product	Scale	Catalog#
LeXHema Panel v2.0, 16 rxn	16 rxn	LX01722
LeXHema Panel v2.0, 96 rxn	96 rxn	LX01721

Statement

For research use only. Not for use in diagnostic procedures.

Without the written permission of LexigenBio, no other individual or organization may reproduce, copy, edit, expand, reduce or translate the contents of this document in any form into other languages for any purpose.

In case of improper use, LexigenBio Inc reserves all rights

LexigenBio Inc. reserves all rights.