

# LeXBot HT




Fully-automated NGS Workstation

# PART 01


## Introduction

**LeXBot HT Fully-automated NGS Workstation**, independently developed by LexigenBio, is specially designed for automated NGS library preparation. It encompasses sample dispensing, nucleic acid extraction, DNA normalization, nucleic acid purification, quantitative loading, and sample mixing, library preparation, and hybrid capture, all aimed at dramatically accelerating laboratory throughput and advancing scientific discovery. The workstation is equipped with a 24-channel pipetting head that flexibly supports the loading of 1/8/16/24 pipette tips and seamlessly integrates a plate gripper, thermal cycler (brand optional), temperature control module, heating and oscillating module, purification module, fluorometer module, UV lamp set, and HEPA filtration system. With 24 built-in SBS standard plate positions, this workstation achieves a maximum throughput of up to 48 pre-libraries or 24 hybrid capture reactions per run, meeting the demands of large-scale, high-throughput applications.



# PART 02

## Application



# PART 03

## Highlight

- 24 standard SBS built-in plate positions, allowing freely configurable dynamic deck layouts to meet diverse experimental needs.
- Seamless integration of multifunctional components, including thermal cycler, temperature control, heating and oscillation, purification, and fluorometer modules.
- Equipped with a HD wide-angle camera for continuous video recording of operations, facilitating retrospective error review.
- A 3D robotic arm with a plate gripper intelligently handles various consumables, ensuring complete deck coverage.

### Flexible Integration



- Multi-channel pipetting (1/8/16/24) with liquid-level detection and intelligent tracking for precise dispensing.
- 24-channel pipetting head offers segmented Z-axis descent (movement, detection, and aspiration).
- Supports a maximum throughput of up to 48 pre-libraries or 24 hybrid capture reactions, meeting varying throughput demands.

### Versatile Application

- Preset for NGS library preparation and hybrid capture; supports custom applications and script development.
- Magnetic sliding window design for effortless operation and enhanced user comfort.
- High-definition touch-sensitive LCD screen with adjustable height and angle for easy viewing and operation.
- Graphical process editing with drag-and-drop script elements enables zero-threshold visual programming for seamless human-machine interaction.

### User-friendly Simplicity




- Bright, multi-color status indicator lights provide real-time updates on operational status, minimizing downtime caused by hardware or software errors.
- Equipped with an emergency pause button for one-touch emergency stops, ensuring human-machine safety.
- Standard HEPA filtration system and UV lamp maintain a continuously clean operating environment.

### Secure and Reliable



# PART 04

## Deck Layout



# PART 05

## Component

### Temperature Control Module



- Consists of 2 modules.
- Temperature control range: 4°C to 75°C.
- Temperature accuracy:  $\leq \pm 1.0^\circ\text{C}$ .
- Temperature uniformity:  $\leq 1.0^\circ\text{C}$ .
- Temperature fluctuation:  $\leq \pm 0.5^\circ\text{C}$ .

### Thermal Cycler Module




- Brand optional.
- Temperature control range: 4°C to 99°C.
- Heated lid temperature range: 30°C to 110°C.
- Average heating rate:  $\geq 3.0^\circ\text{C}/\text{sec}$ .
- Average cooling rate:  $\geq 1.5^\circ\text{C}/\text{sec}$ .
- Temperature accuracy:  $\leq \pm 0.5^\circ\text{C}$ .
- Temperature uniformity:  $\leq 1.0^\circ\text{C}$ .
- Temperature fluctuation:  $\leq \pm 0.2^\circ\text{C}$ .

### Heating and Oscillating Module



- Temperature control range: (room temperature + 5°C) to 99°C.
- Temperature accuracy:  $\leq \pm 2.0^\circ\text{C}$ .
- Temperature uniformity:  $\leq 1.5^\circ\text{C}$ .
- Temperature fluctuation:  $\leq \pm 0.5^\circ\text{C}$ .
- Oscillating speed adjustable from 200 rpm to 2,000 rpm.
- Oscillating speed accuracy:  $\leq \pm 25$  rpm.

### Purification Module



- Function optional.
- movable 96-well PCR plate ring magnet stand.
- liftable 96-well deep-well plate ring magnet stand.

### 24-channel Pipetting Module



- liquid-level detection supporting intelligent liquid-level tracking.
- Segmented Z-axis descent (movement, detection, and aspiration).
- Pipetting range: 2  $\mu\text{L}$  ~ 200  $\mu\text{L}$ .
- Performance:
  - 2  $\mu\text{L}$ : CV < 5%, accuracy <  $\pm 10\%$ .
  - 200  $\mu\text{L}$ : CV < 1%, accuracy <  $\pm 1\%$ .

### Fluorometer Module



- Quantification with support for normalization and pooling.
- Rapid: 3 ~ 5 sec per measurement.
  - Low input: 1  $\mu\text{L}$  ~ 20  $\mu\text{L}$ .
- High sensitivity: 0.01 ng/ $\mu\text{L}$  (dsDNA).
- Detection range: 0.5 ng/ $\mu\text{L}$  ~ 100 ng/ $\mu\text{L}$ .
- Throughput: 8 samples per run, up to 48 samples per batch.

### PCR Plate Carrier



### Composite Plate Carrier



### Reagent Reservoir Rack




### Deep-well Plate Carrier




## PART 06 Consumable


### 200 $\mu\text{L}$ filtered tip



### 50 $\mu\text{L}$ slim-profile filtered tip



### fully skirted 96-well PCR plate



### 1.3 mL U-bottom deep-well plate



### 25 mL reagent reservoir



### 22 mL 12-channel reagent reservoir




### fluorometer 8-strip tube



### 0.2 mL thin-wall PCR 8-strip tube





## PART 07

### Specification

<b>Host Parameters</b>	Size	1,160 mm (width) × 800 mm (depth) × 955 mm (height)		Weight	200 Kg
	Maximum output power	1,700 VA			
<b>Power Requirements</b>	Voltage	AC 100 V ~ 240 V	Frequency	50/60 Hz	
	<b>Operating Environment</b>	Temperature range	10°C ~ 30°C	Relative humidity	20% ~ 80%, non-condensing
Barometric pressure range		86 kPa ~ 106 kPa			
<b>Robot Arm</b>	Repeated positioning	±0.1 mm			

## PART 08

### Ordering Info

Catalog	Product	Scale
LX07504	LeXBot HT Fully-automated NGS Workstation	Unit, RUO

#### Statement

For research use only. Not for use in diagnostic procedures.

Without the written permission of LexigenBio, no other individual or organization may reproduce, copy, edit, expand, reduce or translate the contents of this document in any form into other languages for any purpose.

In case of improper use, LexigenBio Inc reserves all rights

LexigenBio Inc. reserves all rights.



LexigenBio Inc.

✉ [support@lexigenbio.com](mailto:support@lexigenbio.com)  
 🌐 [www.lexigenbio.com](http://www.lexigenbio.com)